

SILICONE HYDROGEL LENSES						
	Colour	Diameter	Base Curve	Power		
O₂OPTIX™ (6-pack) (lotrafilcon B, 33% water)	Handling Tint	14.2	8.6	+6.00 to -10.00 (excludes plano) (0.50D steps above -8.00)		
	Colour	Diameter	Base Curve	Power		
NIGHT & DAY® (6-pack) (lotrafilcon A, 24% water)	Clear	13.8	8.4 8.6	+6.00 to -10.00 (0.50D steps above -8.00)		
FREQUENT REPLACEMENT AND DISPOSABLE LENSES						
Daily Disposable	Colour	Diameter	Base Curve	Cylinder	Power	Axis
Focus® DAILIES® (30 & 90-pk) (nefilcon A, 69% water)	Visitint	13.8	8.6		-0.50 to -10.00 +0.50 to +6.00 (0.50D steps above -6.00)	
Focus DAILIES Progressives (30-pack) (nefilcon A, 69% water)	Visitint	13.8	8.6		+5.00 to -6.00 Progressive add effective up to +3.00D	
Focus DAILIES Toric (30-pk) (nefilcon A, 69% water)	Visitint	14.2	8.6	-0.75	+4.00 to -8.00 (0.50D steps above -6.00)	90°, 180°
FreshLook® ONE-DAY Color Contact Lenses (10-pk) (nefilcon A, 69% water)	Blue, Green, Pure Hazel, Gray	13.8	8.6		Plano -0.50 to -6.00 (0.25D steps)	
Standard Spheres	Colour	Diameter	Base Curve	Power		
Focus Visitint® 1-2 Week (6-pack) (vifilcon A, 55% water)	Visitint	14.0	8.4 8.8	+4.00 to -10.00 (0.50D steps above -6.00)		
Focus Visitint Monthly (6-pack) (vifilcon A, 55% water)	Visitint	14.0	8.6 8.9	+6.00 to -15.00 (0.50D steps above -8.00)		
FreshLook Handling Tint (6-pack) (phemfilcon A, 55% water)	Handling Tint	14.5	Median	+6.00 to -8.00 (0.25D steps) Excludes Plano		
Precision UV® (6-pack) (vasurfilcon A, 74% water)	Handling Tint	14.4	8.4	+8.00 to -10.00 (0.50D steps above +5.00 and -6.00)		
			8.7	+10.00 to -16.00 (0.50D steps above +5.00 and -6.00)		
Cosmetic Spheres	Colour	Diameter	Base Curve	Power		
Focus Softcolors® 1-2 Week (6-pack) (vifilcon A, 55% water)	Royal Blue, Aqua & Evergreen	14.0	8.4 8.8	+4.00 to -6.00 (0.25D steps)		
Focus Softcolors Monthly (6-pack) (vifilcon A, 55% water)	Royal Blue, Aqua & Evergreen	14.0	8.6 8.9	+6.00 to -8.00 (0.25D steps)		
FreshLook Colors (6-pack) (phemfilcon A, 55% water)	Blue, Green, Hazel & Violet	14.5	Median	+6.00 to -8.00 (0.25D steps)		
FreshLook ColorBlends® (6-pack) (phemfilcon A, 55% water)	Blue, Green, Brown, Gray, Honey, Turquoise, Amethyst, Pure Hazel & True Sapphire	14.5	Median	+6.00 to -8.00 (0.25D steps)		
FreshLook Dimensions™ (6-pack) (phemfilcon A, 55% water)	Pacific Blue, Sea Green & Caribbean Aqua	14.5	Median	+6.00 to -8.00 (0.25D steps)		
FreshLook Radiance™ (6-pack) (phemfilcon A, 55% water)	Autumn, Sunrise, Eden & Moonlight	14.5	Median	+6.00 to -8.00 (0.25D steps)		
FreshLook ONE-DAY Color Contact Lenses (10-pk) (nefilcon A, 69% water)	Blue, Green, Pure Hazel, Gray	13.8	8.6		Plano -0.50 to -6.00 (0.25D steps)	
Standard Torics	Colour	Diameter	Base Curve	Cylinder	Power	Axis
Focus Toric Monthly (6-pack) (vifilcon A, 55% water)	Visitint	14.5	8.9	-1.00 & -1.75	plano to -6.00 +0.25 to +4.00	full circle 10° steps 180° ±20° 90° ±20° (in 10° steps)
				-2.50	plano to -6.00	(in 10° steps)
				9.2	-1.00 & -1.75	plano to -6.00
FreshLook Toric (6-pack) (phemfilcon A, 55% water)	Handling Tint	14.5	Median	-0.75, -1.25 & -1.75	plano to -4.00 (0.25D steps)	180° ±20° 90° ±20° (in 10° steps)
Cosmetic Torics	Colour	Diameter	Base Curve	Cylinder	Power	Axis
FreshLook ColorBlends Toric (6-pack) (phemfilcon A, 55% water)	Blue, Green, Honey & Gray	14.5	Median	-0.75, -1.25 & -1.75	plano to -4.00 (0.25D steps)	180° ±20° 90° ±20° (in 10° steps)
Multifocals	Colour	Diameter	Base Curve	Power		
Focus Progressives (6-pack) (vifilcon A, 55% water)	Visitint	14.0	8.6 8.9	+6.00 to -7.00 Single progressive add effective range up to +3.00D		
CONVENTIONAL / SPECIALTY LENSES						
Standard Spheres	Colour	Diameter	Base Curve	MTO+ Power	Std Power	MTO- Power
Cibasoft® Standard Visitint (tefilcon, 37.5% water)	Visitint	13.8	8.3 8.6 8.9		+6.00 to -6.00 (0.25D steps)	
Cibasoft Visitint (tefilcon, 37.5% water)	Visitint	14.5	8.6		plano to -8.00 (0.50D steps above -6.00)	-8.50 to -10.00 (0.50D steps)
			8.9			
			9.2		+6.00 to -8.00 (0.50D steps above -6.00)	-8.50 to -10.00 (0.50D steps)
Clearview® Made-to-Order Only (polymacon, 38% water)	Locator Tint	13.8	8.4	+0.25 to +4.00 (0.25D steps)		plano to -10.00 (0.50D steps above -6.00)
			8.7			
			8.0			
CSI® Daily Wear with Locator Tint (crofilcon A, 38.5% water)	Locator Tint	13.8	8.0	+6.25 to +8.00 (0.25D steps)	+6.00 to -8.00 (0.25D steps)	-8.25 to -20.00 (0.50D steps above -10.00)
			8.3 8.6			

CONVENTIONAL / SPECIALTY LENSES

Standard Spheres Continued						
Standard Spheres Continued	Colour	Diameter	Base Curve	MTO+ Power	Std Power	MTO- Power
CSI Flexible Wear Made-to-Order Only (crofilcon A, 38.5% water)	Clear	14.8	8.6 8.9 9.35			plano to -7.00 (0.25D steps)
		13.8	8.0 8.3 8.6 8.9			plano to -10.00 (0.25D steps)
DuraSoft® 2 LiteTint® (phemfilcon A, 38% water)	Handling Tint	14.5	8.3 8.6 9.0	+6.25 to +20.00 (0.50D steps above +10.00)	+6.00 to -8.00 (0.25D steps)	-8.25 to -20.00 (0.50D steps above -10.00)
		13.8	8.0 8.3 8.6			
DuraSoft 3 LiteTint (phemfilcon A, 55% water)	Handling Tint	14.5	8.3 8.6 9.0	+6.25 to +20.00 (0.50D steps above +10.00)	+6.00 to -8.00 (0.25D steps)	-8.25 to -20.00 (0.50D steps above -10.00)
Hydrocurve® II Made-to-Order Only (bufilcon A, 55% water)	Clear	14.0	8.5	+0.25 to +7.00 (0.50D steps above +4.00)		plano to -12.00 (0.50D steps above -6.00)
		14.5	8.8			
Softmate® II (bufilcon A, 55% water)	Clear	14.3	8.7	+6.50 to +7.00 (0.50D steps)	+6.00 to -8.00 (0.50D steps above +4.00 & -6.00)	-8.50 to -12.00 (0.50D steps)
		14.8	9.0			
Cosmetic Spheres						
Cosmetic Spheres	Colour	Diameter	Base Curve	MTO+ Power	Std Power	MTO- Power
Cibasoft Softcolors (tefilcon, 37.5% water)	Aqua, Blue, Evergreen, Green, Amber, Royal Blue & Violet* ***** * Violet in 13.8 Diameter only	13.8	8.3 8.6 8.9		+6.00 to -8.00 (0.50D steps above -6.00)	-8.50 to -10.00 (0.50D steps)
		14.5	8.6 8.9 9.2		plano to -8.00 (0.50D steps above -6.00)	-8.50 to -10.00 (0.50D steps)
DuraSoft 2 Colors (phemfilcon A, 38% water)	Blue, Green, Hazel & Gray	14.5	8.3	+6.25 to +20.00 (0.50D steps above +10.00)	+6.00 to -8.00 (0.25D steps)	-8.25 to -20.00 (0.50D steps above -10.00)
			8.6	+6.25 to +20.00 (0.50D steps above +10.00)	+6.00 to -8.00 (0.25D steps)	-8.25 to -20.00 (0.50D steps above -10.00)
			9.0	+0.25 to +20.00 (0.50D steps above +10.00)		plano to -20.00 (0.50D steps above -10.00)
DuraSoft 2 Colors for Light Eyes (phemfilcon A, 38% water)	Sky Blue, Aquamarine, Jade Green & Violet-Blue	14.5	8.3	+0.25 to +20.00 (0.50D steps above +10.00)		plano to -20.00 (0.50D steps above -10.00)
			8.6	+6.25 to +20.00 (0.50D steps above +10.00)	+6.00 to -8.00 (0.25D steps)	-8.25 to -20.00 (0.50D steps above -10.00)
			9.0	+0.25 to +20.00 (0.50D steps above +10.00)		plano to -20.00 (0.50D steps above -10.00)
DuraSoft 2 ColorBlends® (phemfilcon A, 38% water)	Blue, Green, Brown & Gray	14.5	8.3 8.6		+4.00 to -6.00 (0.25D steps)	
DuraSoft 3 Colors (phemfilcon A, 55% water)	Baby Blue, Emerald Green, Hazel, Aqua, Sapphire Blue, Jade Green, Violet, Misty Gray, Chestnut Brown & Baby Blue Enhance	14.5	8.3	+6.25 to +20.00 (0.50D steps above +10.00)	+6.00 to -8.00 (0.25D steps)	-8.25 to -20.00 (0.50D steps above -10.00)
			8.6	+6.25 to +20.00 (0.50D steps above +10.00)	+6.00 to -8.00 (0.25D steps)	-8.25 to -20.00 (0.50D steps above -10.00)
			9.0	+0.25 to +20.00 (0.50D steps above +10.00)		plano to -20.00 (0.50D steps above -10.00)
DuraSoft 3 Complements® (phemfilcon A, 55% water)	Blue, Green, Brown, Blue-Violet & Shadow Gray	14.5	8.3	+0.25 to +20.00 (0.50D steps above +10.00)		plano to -20.00 (0.50D steps above -10.00)
			8.6	+6.25 to +20.00 (0.50D steps above +10.00)	+6.00 to -8.00 (0.25D steps)	-8.25 to -20.00 (0.50D steps above -10.00)
			9.0	+0.25 to +20.00 (0.50D steps above +10.00)		plano to -20.00 (0.50D steps above -10.00)
DuraSoft 3 ColorBlends (phemfilcon A, 55% water)	Blue, Green, Brown, Gray, Honey, Turquoise & Amethyst* ***** * Amethyst in 8.6 Base Curve with +6.00 to -8.00 powers only	14.5	8.3	+0.25 to +20.00 (0.50D steps above +10.00)		plano to -20.00 (0.50D steps above -10.00)
			8.6	+6.25 to +20.00 (0.50D steps above +10.00)	+6.00 to -8.00 (0.25D steps)	-8.25 to -20.00 (0.50D steps above -10.00)
			9.0	+0.25 to +20.00 (0.50D steps above +10.00)		plano to -20.00 (0.50D steps above -10.00)
Elegance® (polymacon, 38% water)	Light Blue, Dark Blue, Light Green, Dark Green, Grey, Hazel	13.8	8.4 8.7		+1.00 to -6.00 (0.25D steps)	
Illusions® (tefilcon, 37.5% water)	Violet, Soft Amber, Gray, Deep Blue, Soft Blue, Deep Green & Soft Green	13.8	8.3 8.6 8.9		+4.00 to -8.00 (0.50D steps above -6.00)	-8.50 to -10.00 (0.50D steps)
WildEyes® (phemfilcon A, 55% water)	wildfire, white-out, zoomin' hypnotica, red hot, cat eye, icefire, knockout, black-out, jaguar & zebra	14.5	8.6		plano to -6.00 (0.25D steps)	
Cosmetic Torics						
Cosmetic Torics	Colour	Diameter	Base Curve	Cylinder	Power	Axis
DuraSoft 2 OptiFit Toric Colors Made-to-Order Only (phemfilcon A, 38% water)	Blue, Green, Hazel & Gray	14.5	Steep Median Flat	-0.75 to -9.75 (0.50D steps)	+10.00 to -20.00 (0.50D steps above +4.00 and -8.00)	full circle 5°
DuraSoft 2 OptiFit Toric Colors for Light Eyes Made-to-Order Only (phemfilcon A, 38% water)	Sky Blue, Jade Green, Violet-Blue & Aquamarine	14.5	Steep Median Flat	-0.75 to -9.75 (0.50D steps)	+10.00 to -20.00 (0.50D steps above +4.00 and -8.00)	full circle 5°
DuraSoft 3 OptiFit Toric Colors Made-to-Order Only (phemfilcon A, 55% water)	Baby Blue, Emerald Green, Hazel, Aqua, Sapphire Blue, Jade Green, Violet, Misty Gray, Chestnut Brown & Baby Blue Enhance	14.5	Steep Median Flat	-0.75 to -5.75 (0.50D steps)	+4.00 to -8.00 (0.25D steps)	full circle 5°

CONVENTIONAL / SPECIALTY LENSES

Cosmetic Torics Continued						
Colour	Diameter	Base Curve	Cylinder	Power	Axis	
DuraSoft 3 OptiFit Toric Complements Made-to-Order Only (phemfilcon A, 55% water)	Blue, Green, Brown, Blue-Violet & Shadow Gray	14.5	Steep Median Flat	-0.75 to -5.75 (0.50D steps)	+4.00 to -8.00 (0.25D steps)	full circle 5°
DuraSoft 3 OptiFit ColorBlends Made-to-Order Only (phemfilcon A, 55% water)	Blue, Green, Brown, Gray, Honey, Turquoise & Amethyst	14.5	Steep Median Flat	-0.75 to -5.75 (0.50D steps)	+4.00 to -8.00 (0.25D steps)	full circle 5°
Standard Torics						
Colour	Diameter	Base Curve	Cylinder	Power	Axis	
CSI Toric Made-to-Order Only (crofilcon A, 38.5% water)	Locator Tint	14.0	8.3 8.6	-1.00, -1.75 & -2.50	+ 4.00 to -8.00 (0.25D steps)	full circle 5°
CIBASOFT Progressive Toric Made-to-Order Only (tefilcon, 37.5% water)	Clear	14.5	8.6 8.9	-0.75 to -2.75 (0.25 steps)	+9.00 to -9.00 (0.25 steps)	full circle 5°
Single progressive add effective range up to +3.00D Zero cylinder lens available for spherical patients						
DuraSoft 2 OptiFit Toric Made-to-Order Only (phemfilcon A, 38% water)	Handling Tint	14.5	Steep Median Flat	-0.75 to -9.75 (0.50 steps)	+10.00 to -20.00 (0.50D steps above +4.00 and -8.00)	full circle 5°
DuraSoft 3 OptiFit Toric Made-to-Order Only (phemfilcon A, 55% water)	Handling Tint	14.5	Steep Median Flat	-0.75 to -5.75 (0.50 steps)	+4.00 to -8.00 (0.25D steps)	full circle 5°
DuraSoft 3 OptiFit Toric UV Made-to-Order Only (phemfilcon A, 55% water)	Handling Tint	14.5	Steep Median Flat	-0.75 to -5.75 (0.50D steps)	+4.00 to -8.00 (0.25D steps)	full circle 5°
Hydrocurve 3 Toric Made-to-Order Only (butilcon A, 55% water)	Clear	14.5	8.8	-0.75, -1.25 & -2.00	+4.00 to -8.00 (0.25D steps)	full circle 5°
Torisoft® (tefilcon, 37.5% water)	Clear	14.5	8.9	-1.00 & -1.75	+4.00 to -7.00 (0.50D steps above -6.00)	180° ±30° 90° ±30° (in 10° steps)
				-2.50	plano to -7.00 (0.50D steps above -6.00)	180° ±20° 90° ±20° (in 10° steps)
Torisoft Made-to-Order (tefilcon, 37.5% water)	Clear	14.5	8.9	-1.00 & -1.75	+4.00 to -7.00 (0.50D steps above -6.00)	Oblique axis 40°, 50°, 130° & 140°
			8.6	-1.00 & -1.75	+4.00 to -7.00 (0.50D steps above -6.00)	full circle 10°
			9.2	-1.00 & -1.75	plano to -7.00 (0.50D steps above -6.00)	180° ±20° 90° ±20° (in 10° steps)
Other Specialty Lenses						
Colour	Diameter	Base Curve	MTO+ Power	Std Power	MTO- Power	
Hydrocurve II Aphakic Made-to-Order Only (butilcon A, 55% water)	Clear	14.0	8.5	+7.50 to +20.00 (0.50D steps)		
		14.5	8.8	+7.50 to +20.00 (0.50D steps)		
		15.5	9.2	+12.00 to +16.00 (0.50D steps)		
			9.5 9.8	+7.50 to +20.00 (0.50D steps)		
Hydrocurve II Bifocal Made-to-Order Only (butilcon A, 55% water)	Clear	14.8	9.0	+4.00 to -6.00 (Single progressive add effective up to +1.50D)		
SoftPerm® (RGP center with 25% water hydrophilic skirt)	Clear	14.3	6.5 6.7 6.9 7.1 to 8.1 (in 0.10mm steps)	-3.00 to -16.00 (0.50D steps above -8.00)		
+6.00 to -13.00 (0.50 steps above +4.00 and -8.00)						
PROSTHETIC LENSES						
Prosthetic Lenses	Colour	Diameter	Base Curve	Power		
DuraSoft 2 Prosthetic Made-to-Order (phemfilcon A, 38% water)	Available in all DuraSoft Colors, Complement & ColorBlend colours, Available in 5 print patterns.	13.5, 13.8, 14.0, 14.5, & 15.0	7.80 to 9.40	+20.00 to -20.00 (0.25D Steps) Toric also available; Call for information		
DuraSoft 3 Prosthetic Made-to-Order (phemfilcon A, 55% water)	Available in all DuraSoft Colors, Complement & ColorBlend colours, Available in 5 print patterns.	13.5, 13.8, 14.0, 14.5, & 15.0	7.80 to 9.40	+20.00 to -20.00 (0.25D Steps) Toric also available; Call for information		

Note: For additional information and to order, call the Special Eyes Program consultant at (800) 488-6859.

LENS CARE PRODUCTS

Lens Care Kits	Case/Pack	Contents
Clear Care® Starter Kit	24	120 ml Clear Care Cleaning and Disinfecting Solution 1 disposable Cup and Disc
SOLOCARE AQUA™ All-In-One Solution Starter Kit	24	60 ml SOLOCARE AQUA All-In-One Solution 1 Lens Case (MicroBlock)
Other Lens Care Products for Resale		Case/Pack
DISINFECTANTS		
AOSEPT® PLUS Disinfecting Solution 360 ml	6	
AOSEPT Disposable Cup & Disc	6	
Clear Care 360 ml	6	
Clear Care 240 ml Twin Pack	12	
Clear Care 360 ml Twin Pack	12	
SOLOCARE AQUA All-In-One Solution 360 ml	6	
SOLOCARE AQUA All-In-One Solution 360 ml Twin Pack	12	
SOLOCARE AQUA All-In-One Solution 120 ml	6	
LUBRICANTS & REWETTING DROPS		
AQuify® Long-Lasting Comfort Drops 10 ml	6	
DAILY CLEANERS		
AOFLOW® Extra Strength Daily Cleaner 20 ml	6	
SALINES		
SoftWear® Saline 360 ml	6	
ENZYMATIC CLEANERS		
Unizyme® 12-pack Trade	6	

© 2006 CIBA Vision